

# **Spectral Resolution Add-on for Two Photon Microscope**

**(OTT ID 1105-1242)**

**Inventor: Dr. Valerică Raicu, Physics Department  
UW-Milwaukee**


**For further information please contact:**

**Dr. Smruti Patil  
Licensing Associate  
1440 East North Ave.  
Milwaukee, WI 53202  
Tel: 414-906-4657  
smruti@uwmrf.org**

## Typical two-photon microscopes scramble dynamic spectral information during data acquisition

- ❑ Samples are scanned with a laser in a serial, i.e. point by point, fashion which limits real-time imaging of dynamic processes.
- ❑ Beam splitters and filters, used in these microscopes to separate signal from different fluorescent molecules, have to be changed after each serial scan, and therefore spectral information about dynamic molecular complexes is scrambled before all the spectral information is acquired.
- ❑ Spectral bleed-through, which is the crossing of signal from one fluorescent molecule into the detection channel designated for another, occurs when beam splitters are used to separate signal from different fluorescent molecules.
- ❑ Because of spectral bleed-through, fluorescent molecules have to be chosen that do not spectrally overlap, which severely limits the range of fluorescent molecules which may be employed.

We present a module upgrade which separates fluorescence signal obtained from multiple locations in a sample into hundreds of wavelength channels simultaneously



Fluorescence intensity maps of Chinese Hamster Ovary cells expressing rhodopsin labeled with SCFP3A.

- ❑ Fluorescence emission from a sample is projected through a transmission grating onto an electron-multiplying CCD (EM-CCD) camera with single photon sensitivity. Different wavelength channels are all detected at the same time
- ❑ Because the detector is a 2D array, multiple points in a sample can be excited and imaged in parallel, allowing for a two order of magnitude increase in speed and/or sensitivity

**This module is specifically designed to deliver real-time, spectrally resolved analysis of single molecules and molecular complexes and their spatial distribution in living cells**

# Add-on Device (Specifications)

- Microscope**-Module can be directly coupled to left or right side ports on a wide range of microscopes sold by any of the major manufacturers.
- Detector**-Either the Andor iXon X3 or iXon Ultra 897 EMCCD can used for detection of emitted fluorescence
- Laser**-Can be coupled to any femtosecond pulsed free-space laser emitting above 680 nm
- Computer**-Compatible with any standard PC running either Windows 7 or Windows 10. A National Instruments PCIE DAQ Card must be installed on the computer



## Applications

- Live Cell Imaging
- 3D Imaging
- Deep Tissue Imaging
- Long Term Imaging without compromised tissue viability
- Whole Organ or Slice Imaging
- Many more

## US and European issued patents

US7,973,927, US8,094,304, US9,103,721, US8,982,206 and EP20120782796 validated in Belgium, France, Germany, United Kingdom, and Ireland

## Current Status

- This technology can be sold separately or as an add-on attachment with new microscope
- Prof. Raicu's lab also offers this technology on service basis along with operating systems
- Technology is available for licensing under exclusive or non-exclusive terms

- ❑ Characterizing the effect of ligands on protein-protein interactions, particularly on association of membrane receptors using 2D FIF method.
  - The method is based on 2D FIF, published by Stoneman et al in Nature Methods in 2019
  
- ❑ Screening of chemical compounds based on their ability to bind to the drug target (i.e., cell cultures expressing receptors of interest)
  
- ❑ Cost of service will depend on type of sample and sample size, and is based on an hourly fee for the imaging facility and technician/scientist time

- ❑ The Raicu lab has developed a two-photon microscope add-on device that increases both resolution and data acquisition speed.
- ❑ The add-on design is complete and several have been manufactured. Operating software has also been written and is available.
- ❑ The technology has been used in a number of quantitative FRET imaging publications where the binding stoichiometry and geometry of protein complexes has been determined

# **Spectral Resolution Add-on for Two Photon Microscope**

**(OTT ID 1105-1242)**

**For further information please contact:**

**Dr. Smruti Patil**  
**Licensing Associate**  
**1440 East North Ave.**  
**Milwaukee, WI 53202**  
**Tel: 414-906-4657**  
**[smruti@uwmmrf.org](mailto:smruti@uwmmrf.org)**